



7.972 [202.5] SLIT WIDTH

7.776 [197.5] ARTWORK COLOR WIDTH

7.500 [190.5] ARTWORK COPY WIDTH

4.547 [115.5] CUT LENGTH

4.390 [111.5] COLOR HEIGHT

4.272 [108.5] COPY HEIGHT

IMMUNE SUPPORT • BRAIN DEVELOPMENT • DIGESTIVE HEALTH

Abbott

NEW!

THE PERFECT BOTTLE®
Simply attach nipple and feed

Similac®
360 Total Care®

5 HMO PREBIOTICS

Our closest formula to breast milk

HMOs ARE IMMUNE-NOURISHING PREBIOTICS
STRUCTURALLY IDENTICAL TO THOSE IN BREAST MILK

DHA, LUTEIN, VITAMIN E • NO ARTIFICIAL GROWTH HORMONES*
NON GMO† • NO PALM OLEIN OIL

8 FL OZ (237 mL)

MILK-BASED READY TO FEED

DO NOT ADD WATER

INFANT FORMULA WITH IRON

NOT FOR INFANTS OR CHILDREN WITH GALACTOSEMIA

NUTRIENTS per 100 Calories (5 fl oz)

PROTEIN	2.07 g	WATER	133 g
FAT	5.40 g	LINOLEIC ACID	1000 mg
CARBOHYDRATE	11.0 g		

VITAMINS

VITAMIN A	300 IU	NIACIN	1050 mcg
VITAMIN D	60 IU	FOLIC ACID (FOLACIN)	15 mcg
VITAMIN E	1.5 IU	PANTOTHENIC ACID	450 mcg
VITAMIN K	8 mcg	BIOTIN	4.4 mcg
THIAMIN (MT. B ₁)	100 mcg	VITAMIN C (ASCORBIC ACID)	9 mg
RIBOFLAVIN (MT. B ₂)	150 mcg	CHOLINE	24 mg
VITAMIN B ₆	60 mcg	INOSITOL	24 mg
VITAMIN B ₁₂	0.25 mcg		

MINERALS

CALCIUM	78 mg	COPPER	90 mcg
PHOSPHORUS	42 mg	IODINE	15 mcg
MAGNESIUM	6 mg	SELENIUM	2 mcg
IRON	1.8 mg	SODIUM	24 mg
ZINC	0.75 mg	POTASSIUM	105 mg
MANGANESE	5 mcg	CHLORIDE	65 mg

INGREDIENTS: WATER, NONFAT MILK, LACTOSE, HIGH OLEIC SAFFLOWER OIL, SOY OIL, COCONUT OIL, WHEY PROTEIN CONCENTRATE, LESS THAN 0.5% OF: HUMAN MILK OLIGOSACCHARIDES*, 2'-FUCOSYLACTOSE, LACTO-N-ETRAOSE, 3-FUCOSYLACTOSE, 6'-SIALYLACTOSE, 3'-SIALYLACTOSE, C. COHNII OIL, M. ALPINA OIL, BETA-CAROTENE, LUTEN, ASCORBIC ACID, SOY LECITHIN, MONOGLYCERIDES, POTASSIUM CITRATE, CALCIUM CARBONATE, POTASSIUM CHLORIDE, CARRAGEENAN, FERROUS SULFATE, MAGNESIUM CHLORIDE, CHOLINE CHLORIDE, CHOLINE BITARTRATE, TAURINE, INOSITOL, CALCIUM PHOSPHATE, ZINC SULFATE, POTASSIUM PHOSPHATE, D-ALPHA-TOCOPHEROL, ACETATE, NIACINAMIDE, CALCIUM PANTOTHENATE, L-CARNITINE, VITAMIN A PALMITATE, COPPER SULFATE, THIAMINE HYDROCHLORIDE, RIBOFLAVIN, PYRIDOXINE HYDROCHLORIDE, FOLIC ACID, MANGANESE SULFATE, PHYLLOQUINONE, BIOTIN, SODIUM SELENATE, VITAMIN D₃, VITAMIN B₁₂, SALT, POTASSIUM IODIDE, POTASSIUM HYDROXIDE, AND NUCLEOTIDES (ADENOSINE 5'-MONOPHOSPHATE, CYTIDINE 5'-MONOPHOSPHATE, DISODIUM GUANOSINE 5'-MONOPHOSPHATE, DISODIUM URIDINE 5'-MONOPHOSPHATE), CONTAINS MILK AND SOY INGREDIENTS.

Abbott Nutrition, Abbott Laboratories
Columbus, Ohio 43219-3034 USA

HMO, NOT FROM HUMAN MILK. SOURCE OF DHA, IS SOURCE OF ARA.

GLUTEN FREE

*NO SIGNIFICANT DIFFERENCE HAS BEEN SHOWN BETWEEN MILK DERIVED FROM rHST-TREATED AND NON-rHST-TREATED COWS

†INGREDIENTS NOT GENETICALLY ENGINEERED

USE BY DATE ON BOTTLE

0 170074 68140 5

USE AS DIRECTED BY A DOCTOR
Directions for Preparation and Use

DO NOT ADD WATER • SHAKE VERY WELL
DO NOT USE if breakaway ring is missing or broken

Your baby's health depends on carefully following these directions. Ask your baby's doctor if you need to boil (sterilize) nipples and rings before use.

1 TWIST OFF cap

2 TWIST ON clean nipple and ring

For use in a different feeding container (cup or bottle):

- Ask your baby's doctor if you need to boil (sterilize) the empty bottles before adding formula.
- Pour desired amount into clean feeding container. Recap bottle, refrigerate and use within 48 hours.

Once feeding begins, use within 1 hour or discard.

Storage: Store unopened at room temperature; avoid extreme temperatures. Do not reuse this bottle. Protect from light.

Warning: Never use a microwave to warm formula. Serious burns can result.

THIS UNIT NOT INTENDED FOR INDIVIDUAL SALE

Pat. www.abbott-usa.com/patents © 2020 Abbott Laboratories 68138_C6655-01

100mm

90mm

85mm

77.5mm

100% 2.8937 [73.5]

4mm

0mm

135% 0.4921 [12.5]

119% 0.3937 [10.0]

100% 0.1969 [5.0]

111% 0.2953 [7.5]

742 DISTORTIONS WILL OCCUR IN AREAS WHERE THERE IS A PERCENTAGE MARKED. PULL OUT GRAPHICS AND COPY AT THE PERCENTAGE LISTED BELOW. HORIZONTALLY ONLY.

Abbott Nutrition

68138-8-SL-C6655-3TC-SWA-D

PROCESS CYAN

PROCESS MAGENTA

PROCESS YELLOW

PROCESS BLACK

PMS 287 C

WHITE

SC SLIP COAT

V VARNISH

PROOF:

EPSON

KODAK

APPROVED*

DISAPPROVED

PRINTER: AMERICAN FUJI SEAL

CAD: PK103212 REV G
PREV. COMMODITY: N/A

GS†: GL0330
†Generic Specification

SIGNATURE

*Approved proof is accurate for copy, size, upc, and process builds. Match spot color(s) to approved ink drawdown(s) if applicable.